

BIOXIS

PHARMACEUTICALS

#BIOXIS BEYOND
INNOVATION



MEDIUM

Is indicated for:

The filling of moderate cutaneous depressions as perioral or commissural area

HA concentration: 25 mg/ml

Volume: 1.1ml


Needle size: 27^{1/2}G or 30G

Injection: Mid and Upper dermis, lips

Ejection force : 10 - 25 N

BDDE rate : 1.15%

BDDE residual rate : *idem* BDDE residual



CYTOSIAL MEDIUM
is indicated for:

Volume restoration and the filling of cutaneous depressions and wrinkles. Facial reconstruction of volumes in the event of light-to-moderate lipodystrophies, i.e. due to HIV treatment. Cytosial Medium is indicated for the filling of moderate cutaneous depressions as perioral or commissural area.

PERIORAL WRINKLES

LIP VOLUME AND CONTOURS



MANUFACTURING PROCESSES
CE MARK CERTIFIED AND
UNDER CFDA INSPECTION



CERTIFICATION ISO 13485
AND CE MARK BY BSI, THE
HIGHEST EUROPEAN STANDARDS



DO YOU WANT TO KNOW WHAT
BIOXIS PHARMACEUTICALS CAN OFFER YOU?
CONTACT OUR EXCLUSIVE DISTRIBUTOR
TO FIND THE MOST ADAPTED SOLUTION.

contact@bioxis.com

BIOXIS PHARMACEUTICALS

317 avenue Jean Jaures,
69007 Lyon, France
www.bioxis.com

CYTOSIAL®

*CYTOSIAL is a registered trademark
of BIOXIS Pharmaceuticals S.A.S.
CYTOSIAL is manufactured
by BIOXIS Pharmaceuticals S.A.S.*

Your exclusive Distributor:

