

#BIOXIS BEYOND INNOVATION

DEEP WRINKLES

Correct deep lines:

Deep wrinkles Deep nasolabial fold Deep marionette fold

HA concentration: 25 mg/ml

Volume: 1.1ml


Needle size: 271/2G

Injection: Deep and mid dermis

Ejection force: 10-25 N

BDDE rate : 1.50%

BDDE residual rate: < 2 ppm







CYTOSIAL DEEP WRINKLES is indicated for:

Volume restoration and the filling of cutaneous depressions and wrinkles Facial reconstruction of volumes in the event of light-to-moderate lipodystrophies, i.e dueto HIV treatment.

Cytosial Deep Wrinkles is indicated for the filling of deep facial wrinkles, such as nasolabial foldsand marionette lines.

marionette lines



MANUFACTURING PROCESSES CE MARK CERTIFIED AND UNDER CFDA INSPECTION



CERTIFICATION ISO 13485 AND CE MARK BY BSI, THE <u>HIGHEST EUROPEAN</u> STANDARDS



DO YOU WANT TO KNOW WHAT BIOXIS PHARMACEUTICALS CAN OFFER YOU? CONTACT OUR EXCLUSIVE DISTRIBUTOR TO FIND THE MOST ADAPTED SOLUTION.

contact@bioxis.com

BIOXIS PHARMACEUTICALS

317 avenue Jean Jaures, 69007 Lyon, France

www.bioxis.com

CYTOSIAL

CYTOSIAL is a registered trademark of BIOXIS Pharmaceuticals S.A.S. CYTOSIAL is manufactured by BIOXIS Pharmaceuticals S.A.S. Your exclusive Distributor:

